

# Chromatography Columns

**SGL - Scientific Glass Laboratories Ltd** of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

*Providing Scientific Glassware to Laboratories for over 40 years*

**Scientific Glass  
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

[www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Chromatography Columns



## Chromatography Columns with Sintered Disc

Item	Sinter Dia	Socket Size	Length mm	Price Each
R1/10	10	14/23	100	22.50
R1/20	10	14/23	200	22.50
R1/30	10	14/23	300	22.50
R1/40	10	14/23	400	22.50
R2/10	20	19/26	100	23.20
R2/20	20	19/26	200	24.00
R2/30	20	19/26	300	24.00
R2/40	20	19/26	400	24.00
R2/50	20	19/26	500	24.00
R3/30	25	24/29	300	24.55
R3/40	25	24/29	400	24.55
R3/50	25	24/29	500	24.55
R5/40	30	34/35	400	28.30
R5/50	30	34/35	500	28.30
R6/40	40	40/45	400	38.40
R6/50	40	40/45	500	38.40

Columns 10mm to 25mm bore fitted with 2mm bore stopcocks.

Columns 30mm to 40mm bore fitted with 4mm bore stopcocks.

Glass or PTFE s/cs - please specify - /G or /P

Por 0, 1, 2, 3 and 4 available.

Add to code no. for example R1/10/0 etc.

E&OE Prices Subject to Change Without Notice May 2018.  
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)

