

Condensers

SGL - Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

Providing Scientific Glassware to Laboratories for over 40 years

**Scientific Glass
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Condensers



Condensers, Air

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CA12/00	30	100	120	B10	B10	6.36
CA12/11	100	200	225	B14	B14	6.02
CA12/12	100	150	170	B19	B19	6.65
CA12/22	200	300	320	B19	B19	7.03
CA12/23	400	400	430	B19	B24	7.00
CA12/33	400	550	570	B24	B24	10.63
CA12/44	500	550	570	B29	B29	11.26

Condensers, Allihn

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CAL20/2	130	200	285	-	B19	25.95
CAL20/3	130	200	285	-	B24	27.23
CAL20/4	130	200	295	-	B29	27.63
CAL20/5	130	200	310	-	B34	25.68
CAL20/6	130	200	310	-	B40	29.77
CAL20/8	150	200	310	-	B50	32.31
CAL30/33	180	30	400	B24	B24	34.55

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.

Condensers, Coil (Graham)

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CC15/22	230	150	265	B19	B19	19.31
CC15/33	230	150	265	B24	B24	20.37
CC30/22	450	300	420	B19	B19	33.42
CC30/33	450	300	420	B24	B24	33.42
CC30/44	450	300	420	B29	B29	33.42
CCX15/4	230	150	260	-	B29	20.80
CCX15/5	230	150	260	-	B34	21.77
CCX20/6	350	200	300	-	B40	28.46
CCX20/7	400	200	315	-	B45	32.51
CCX20/8	400	200	315	-	B50	32.51

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.

E&OE Prices Subject to Change Without Notice May 2018.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Condensers



Condensers, Double Surface (Davies)

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CDS15/11	180	150	240	B14	B14	24.92
CDS15/22	180	150	250	B19	B19	24.67
CDS15/33	230	150	250	B24	B24	27.65
CDS20/22	250	200	300	B19	B19	26.07
CDS20/33	300	200	300	B24	B24	27.97
CDS20/44	300	200	315	B29	B29	30.27
CDS20/3	300	200	290	-	B24	27.02

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.



Condensers, Jacketed Coil

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CCJ15/22	430	150	320	B19	B19	48.48
CCJ15/33	430	150	320	B24	B24	45.53
CCJ15/44	430	150	320	B29	B29	48.50
CCJ15/05	430	150	300	-	B34	45.53

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.



Condensers, Liebig

Item	Effective Length mm	Socket Size	Cone Size	Price Each
CL60/00	60	B10	B10	14.87
CL15/11	150	B14	B14	10.61
CL20/22	200	B19	B19	11.49
CL20/33	200	B24	B24	13.75
CL40/22	400	B19	B19	12.29
CL40/33	400	B24	B24	14.12
CL40/44	400	B29	B29	23.32

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.

E&OE Prices Subject to Change Without Notice May 2018.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Scientific Glass Laboratories
- Traditional Craftsmanship
- Proven Quality

Product Literature available in this series:

- Adapters
- Beakers & Flasks
- Bends
- Chromatography Columns
- Cold Traps
- Condensers
- Dean & Stark
- Digestion Tubes
- Distillation
- Drechsel Bottles & Heads
- Extractors
- Flasks
- Funnels
- Funnels Dropping Separating
- Galileo Liquid Thermometer
- Pipettes, Flasks & Cylinders
- Screw Caps, Stoppers & Stopcocks
- Screwthreads, Discs & Joints
- Specimen Tubes & Vials
- Stirrers
- Tubes & Thermometer Pockets

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk

