

# Thermometers & Tubes

**SGL - Scientific Glass Laboratories Ltd** of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

*Providing Scientific Glassware to Laboratories for over 40 years*

**Scientific Glass  
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

[www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Tubes



## Test Tubes

Item	Capacity ml	Height mm	Dia mm	Socket Size	Price Each
TT6/0	2ml	65	14	B10	3.41
TT10/0	5ml	100	14	B10	3.62
TT10/1	11ml	100	18	B14	3.01
TT10/2	19ml	100	23	B19	3.76
TT12/1	15ml	125	18	B14	3.24
TT15/1	19ml	150	18	B14	2.90
TT15/2	31ml	150	23	B19	3.66
TT15/3	54ml	150	29	B24	3.50
TT20/3	75ml	200	29	B24	3.98
TT25/5	195ml	250	41	B34	9.18

## Tubes, Drying, Angled

Item	Cone Size	Price Each
DTS/0	B10	7.05
DTS/1	B14	7.05
DTS/2	B19	7.93
DTS/3	B24	9.38
DTS/4	B29	9.93

## Tubes / Air Inlets, Bent with Cone

Item	Stem Length	Cone Size	Price Each
MCS/1B	150mm	B14	6.04
MCS/2B	150mm	B19	6.66
MCS/3B	150mm	B24	7.55

## Tubes / Air Inlets, Straight with Cone

Item	Stem Length	Cone Size	Price Each
MCS/0	100mm	B10	5.17
MCS/1	100mm	B14	5.17
MCS/2	100mm	B19	5.42
MCS/3	100mm	B24	5.47
MCS/4	100mm	B29	6.42

E&OE Prices Subject to Change Without Notice May 2018.  
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Thermometer



## Thermometer Pockets

Item	Stem Length	Cone Size	Price Each
TPS/1	45mm	B14	5.40
TPS/2	200mm	B19	5.42

## Thermometers with Cone

Item	Temp Range	Cone Size	Price Each
T1/1	0 -110°C	B14	31.09
T1/2	0 -200°C	B14	31.09
T1/3	0 -300°C	B14	27.95
T1/4	0 -400°C	B14	29.35

E&OE Prices Subject to Change Without Notice April 2019.  
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)

