

Chromatography Columns

SGL - Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

Providing Scientific Glassware to Laboratories for over 40 years

**Scientific Glass
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Chromatography Columns



Chromatography Columns with Sintered Disc

| Item | Sinter Dia | Socket Size | Length mm | Price Each |
|-------|------------|-------------|-----------|------------|
| R1/10 | 10 | 14/23 | 100 | 22.22 |
| R1/20 | 10 | 14/23 | 200 | 22.22 |
| R1/30 | 10 | 14/23 | 300 | 22.22 |
| R1/40 | 10 | 14/23 | 400 | 22.22 |
| R2/10 | 20 | 19/26 | 100 | 22.77 |
| R2/20 | 20 | 19/26 | 200 | 23.57 |
| R2/30 | 20 | 19/26 | 300 | 23.57 |
| R2/40 | 20 | 19/26 | 400 | 23.57 |
| R2/50 | 20 | 19/26 | 500 | 23.57 |
| R3/30 | 25 | 24/29 | 300 | 24.12 |
| R3/40 | 25 | 24/29 | 400 | 24.12 |
| R3/50 | 25 | 24/29 | 500 | 24.12 |
| R5/40 | 30 | 34/35 | 400 | 27.81 |
| R5/50 | 30 | 34/35 | 500 | 27.81 |
| R6/40 | 40 | 40/45 | 400 | 37.71 |
| R6/50 | 40 | 40/45 | 500 | 37.71 |

Columns 10mm to 25mm bore fitted with 2mm bore stopcocks.

Columns 30mm to 40mm bore fitted with 4mm bore stopcocks.

Glass or PTFE s/cs - please specify - /G or /P

Por 0, 1, 2, 3 and 4 available.

Add to code no. for example R1/10/0 etc.

E&OE Prices Subject to Change Without Notice May 2017.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk

