

Digestion Tubes

SGL - Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

Providing Scientific Glassware to Laboratories for over 40 years

**Scientific Glass
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Digestion Tubes



Digestion Tubes, Plain

Item	Capacity ml	Dimensions mm	Price Each
DT1/75	75ml	255 x 26	4.03
DT1/100	100ml	300 x 26	4.13
DT1/250	250ml	300 x 42	6.39

Digestion Tubes, with Joint

Item	Capacity ml	Dimensions mm	Socket Size	Price Each
DT2/15	15ml	145 x 16	B14	4.03
DT2/75	75ml	255 x 26	B24	4.73
DT2/100	100ml	300 x 26	B24	4.81
DT2/250	250ml	300 x 42	B40	6.94

Digestion Tubes, with Spherical Joint

Item	Capacity ml	Dimensions mm	Price Each
DT3/75	75ml	250 x 26	8.76
DT3/100	100ml	300 x 26	8.98
DT3/100/G	100ml	300 x 26 graduated	13.56
DT3/250	250ml	280 x 42	12.35

Digestion Tubes, with Constriction, Ungraduated

Item	Capacity ml	Dimensions mm	Price Each
DT4/75	75ml	255 x 26	6.28
DT4/100	100ml	300 x 26	7.46
DT4/250	250ml	300 x 42	10.50

E&OE Prices Subject to Change Without Notice May 2017.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk

