

Adapters

SGL - Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

Providing Scientific Glassware to Laboratories for over 40 years

**Scientific Glass
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Adapters



Adapters, Cone / Rubber Tubing, right Angled Connection

Item	Cone Size	Price Each
ATA/0	B10	4.74
ATA/1	B14	4.38
ATA/2	B19	4.60
ATA/3	B24	4.64
ATA/4	B29	5.72
ATA/5	B34	6.68

For GL14 thread type add /T to code no.



Adapters, Cone / Rubber Tubing, Straight Connection

Item	Cone Size	Price Each
ATS/1	B14	4.38
ATS/2	B19	4.60
ATS/3	B24	4.64

For GL14 thread type add /T to code no.



Expansion Adapters

Item	Cone Size	Socket Size	Price Each
AX/10	B10	B14	5.48
AX/21	B14	B19	4.34
AX/31	B14	B24	3.76
AX/32	B19	B24	3.76
AX/41	B14	B29	5.80
AX/42	B19	B29	5.80
AX/43	B24	B29	4.78
AX/52	B19	B34	6.74
AX/53	B24	B34	5.74
AX/54	B29	B34	6.22
AX/64	B29	B40	8.48

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Multiple, Two Necks Parallel

Item	Socket Size	Cone Size	Price Each
AMP/00	B10	B10	9.90
AMP/11	B14	B14	9.52
AMP/22	B19	B19	9.90
AMP/32	B19	B24	10.34
AMP/33	B24	B24	12.98
AMP/44	B29	B29	15.02
AMP/52	B19	B34	10.34

Adapters, Multiple, One Parallel & One 45 Degree Neck

Item	Socket Size	Cone Size	Price Each
AMA/22	B19	B19	10.18
AMA/32	B19	B24	10.18
AMA/33	B24	B24	10.18
AMA/43	B24	B29	10.64

Adapters, Multiple, Two Parallel & One 45 Degree Neck

Item	Socket Size	Cone Size	Price Each
AM/322	B19	B24	15.02
AM/422	B19	B29	15.02

Adapters, Receiver, Intermediate (Perkin) All Glass Taps

Item	Socket Size	Cone Size	Price Each
RIP/11	B14	B14	107.44
RIP/22	B19	B19	112.98
RIP/23	B19	B24	112.98
RIP/33	B24	B24	112.06

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Vertical, Sloping

Item	Socket Size	Cone Size	Price Each
ARA/22B	B19	B19	9.90

Adapters, Vertical, Delivery

Item	Socket Size	Cone Size	Price Each
ARA/22	B19	B19	7.08

Adapters, Receiver, Vacuum (to take 4 x 100ml flasks)

Item	Socket Size	Cone Size	Price Each
ARA/21	B19	B14	24.82

Adapters, Receiver, Vertical with Vacuum Connection

Item	Socket Size	Cone Size	Price Each
RASV/22	B19	B19	10.68
RASV/33	B24	B24	10.68

For GL14 thread type add /T to code no.

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Reduction Adapters

Item	Socket Size	Cone Size	Price Each
AR/01	B10	B14	3.36
AR/12	B14	B19	2.92
AR/13	B14	B24	3.20
AR/14	B14	B29	3.46
AR/23	B19	B24	3.04
AR/24	B19	B29	3.42
AR/25	B19	B34	3.64
AR/26	B19	B40	5.56
AR/27	B19	B45	6.16
AR/34	B24	B29	3.36
AR/35	B24	B34	3.52
AR/36	B24	B40	5.20
AR/37	B24	B45	5.46
AR/38	B24	B50	6.68
AR/39	B24	B55	6.90
AR/45	B29	B34	4.22
AR/47	B29	B45	6.46
AR/56	B34	B40	5.68
AR/58	B34	B50	7.20

Adapters, Cone / Screw Thread

Item	Cone Size	Thread Size	Dia Through Hole mm	Price Each
AST/1	B14	14	6.0 - 7.0	4.64
AST/2	B19	14	6.0 - 7.0	4.64
AST/3	B24	14	6.0 - 7.0	5.80
AST/4	B29	14	6.0 - 7.0	6.62
AST/5	B34	14	6.0 - 7.0	7.44
AST/6	B14	18	7.0 - 8.5	5.00
AST/7	B19	18	7.0 - 8.5	4.64
AST/8	B24	18	7.0 - 8.5	6.62
AST/9	B29	18	7.0 - 8.5	6.62
AST/11	B19	24	10.5 - 11.5	7.18
AST/12	B24	24	10.5 - 11.5	6.62
AST/13	B29	24	10.5 - 11.5	7.18
AST/14	B34	24	10.5 - 11.5	8.88

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Cone / Socket with 'T' Connection

Item	Socket Size	Cone Size	Price Each
ATR/11	B14	B14	5.20
ATR/12	B14	B19	5.38
ATR/13	B14	B24	5.80
ATR/22	B19	B19	7.18
ATR/23	B19	B24	5.40
ATR/24	B19	B29	6.06
ATR/25	B19	B34	6.26
ATR/33	B24	B24	7.20
ATR/44	B29	B29	7.64

For GL14 thread type add /T to code no.

Adapters, Socket / Rubber Tubing, Straight Connection

Item	Socket Size	Price Each
ATF/2	B19	5.64
ATF/3	B24	5.64
ATF/4	B29	5.80

For GL14 thread type add / T to code no.

Adapters, Swan Neck, Thread Size GL14

Item	Socket Size	Cone Size	Price Each
ASN/1	B14	B14	9.62
ASN/2	B19	B19	9.20
ASN/3	B19	B24	9.36
ASN/4	B24	B24	11.34

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Stopcock, Cone / Rubber Tubing, Straight Connection

Item	Cone Size	Price Each
ASS/1	B14	14.08
ASS/2	B19	13.26
ASS/3	B24	13.14
ASS/4	B29	16.32

For GL14 thread type add /T to code no.
For P.T.F.E key add P to code no.

Adapters, Stopcock, Cone / Rubber Tubing, Bent Connection

Item	Cone Size	Price Each
ASA/1	B14	15.90
ASA/2	B19	15.90
ASA/3	B24	16.06
ASA/4	B29	16.36

For GL14 thread type add /T to code no.
For P.T.F.E key add P to code no.

Adapters, Stopcock, Socket / Rubber Tubing, Straight Connection

Item	Cone Size	Price Each
ASF/1	B14	15.90
ASF/2	B19	15.90
ASF/3	B24	16.36
ASF/2B	B19	16.06

For GL14 thread type add /T to code no.
For P.T.F.E key add P to code no.

Adapters, Straight Delivery

Item	Socket Size	Price Each
RAS/1	B14	4.44
RAS/2	B19	4.66
RAS/3	B24	4.90
RAS/4	B29	5.10
RAS/5	B34	5.28

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk

