

Condensers

SGL - Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

Providing Scientific Glassware to Laboratories for over 40 years

**Scientific Glass
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Condensers



Condensers, Air

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CA12/00	30	100	120	B10	B10	4.54
CA12/11	100	200	225	B14	B14	4.26
CA12/12	100	150	170	B19	B19	4.80
CA12/22	200	300	320	B19	B19	5.12
CA12/23	400	400	430	B19	B24	5.10
CA12/33	400	550	570	B24	B24	8.24
CA12/44	500	550	570	B29	B29	8.78



Condensers, Allihn

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CAL20/2	130	200	285	-	B19	20.18
CAL20/3	130	200	285	-	B24	22.02
CAL20/4	130	200	295	-	B29	22.34
CAL20/5	130	200	310	-	B34	20.76
CAL20/6	130	200	310	-	B40	24.08
CAL20/8	150	200	310	-	B50	26.12
CAL30/33	180	30	400	B24	B24	27.94

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.

Condensers, Coil (Graham)

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CC15/22	230	150	265	B19	B19	14.94
CC15/33	230	150	265	B24	B24	15.76
CC30/22	450	300	420	B19	B19	26.36
CC30/33	450	300	420	B24	B24	26.08
CC30/44	450	300	420	B29	B29	15.14
CCX15/4	230	150	260	-	B29	16.10
CCX15/5	230	150	260	-	B34	16.84
CCX20/6	350	200	300	-	B40	22.02
CCX20/7	400	200	315	-	B45	25.14
CCX20/8	400	200	315	-	B50	25.14

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.



Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Condensers



Condensers, Double Surface (Davies)

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CDS15/11	180	150	240	B14	B14	19.96
CDS15/22	180	150	250	B19	B19	19.76
CDS15/33	230	150	250	B24	B24	22.14
CDS20/22	250	200	300	B19	B19	20.88
CDS20/33	300	200	300	B24	B24	22.40
CDS20/44	300	200	315	B29	B29	24.24
CDS20/3	300	200	290	-	B24	21.64

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.



Condensers, Jacketed Coil

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CCJ15/22	430	150	320	B19	B19	38.82
CCJ15/33	430	150	320	B24	B24	36.46
CCJ15/44	430	150	320	B29	B29	38.84
CCJ15/05	430	150	300	-	B34	36.46

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.



Condensers, Liebig

Item	Effective Length mm	Socket Size	Cone Size	Price Each
CL60/00	60	B10	B10	12.60
CL15/11	150	B14	B14	8.99
CL20/22	200	B19	B19	9.73
CL20/33	200	B24	B24	11.66
CL40/22	400	B19	B19	10.42
CL40/33	400	B24	B24	11.97
CL40/44	400	B29	B29	19.76

GL14 threads instead of W/connections available at no extra cost. Add /T to code no.

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk

