

Flasks

SGL - Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

Providing Scientific Glassware to Laboratories for over 40 years

**Scientific Glass
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Flasks



Flasks, Buchner, Filter

Item	Capacity ml	Socket Size	Height mm	Dia mm	Price Each
FCB100/B19	100ml	19/26	105	70	9.20
FCB100/B24	100ml	24/29	105	70	9.64
FCB250/B24	250ml	24/29	140	90	10.64
FCB250/B29	250ml	29/32	140	90	11.76
FCB500/B24	500ml	24/29	175	110	12.40
FCB500/B29	500ml	29/32	175	110	12.40
FCB1L/B24	1000ml	24/29	220	135	16.92
FCB1L/B29	1000ml	29/32	220	135	17.72

For GL thread type add /T to code no.



Flasks, Conical Erlenmeyer, Standard Wall

Item	Capacity ml	Socket Size	Price Each
FC10/B14	10ml	B14	4.88
FC25/B10	25ml	B10	4.58
FC25/B14	25ml	B14	4.60
FC25/B19	25ml	B19	4.94
FC50/B14	50ml	B14	4.58
FC50/B19	50ml	B19	4.42
FC50/B24	50ml	B24	3.88
FC100/B14	100ml	B14	5.08
FC100/B19	100ml	B19	3.92
FC100/B24	100ml	B24	4.02
FC100/B29	100ml	B29	5.46
FC150/B19	150ml	B19	5.70
FC150/B24	150ml	B24	4.70
FC150/B29	150ml	B29	5.90
FC250/B19	250ml	B19	4.74
FC250/B24	250ml	B24	4.30
FC250/B29	250ml	B29	5.22
FC250/B34	250ml	B34	5.48
FC500/B19	500ml	B19	5.70
FC500/B24	500ml	B24	4.80
FC500/B29	500ml	B29	5.90
FC500/B34	500ml	B34	6.02
FC1L/B24	1000ml	B24	6.14
FC1L/B29	1000ml	B29	7.60
FC1L/B34	1000ml	B34	7.90
FC2L/B29	2000ml	B29	12.20
FC2L/B34	2000ml	B34	13.70

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Flasks



Iodine Flasks

Item	Capacity ml	Socket Size	Price Each
FI250/B24	250ml	B24	12.62
FI250/B29	250ml	B29	13.90
FI500/B24	500ml	B24	14.30
FI500/B29	500ml	B29	15.10

Iodine Stopper

Item	Cone Size	Price Each
FIS/B24	B24	5.74
FIS/B29	B29	5.74

Kjeldahl Flasks

Item	Capacity ml	Socket Size	Price Each
FK100/B24	100ml	B24	6.06
FK300/B24	300ml	B24	6.36
FK500/B24	500ml	B24	6.36
FK800/B29	800ml	B24	8.80

Round Bottom Flasks (Kjeldahl)

Item	Capacity ml	Socket Size	Price Each
FMC22/50	50ml	B19	4.64
FMC22 100	100ml	B19	5.10

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Flasks



Pear Shaped Flasks, Claisen

Item	Capacity ml	Cone Size	Socket Size	Price Each
FPC3/10	10ml	B14	B14	10.44
FPC3/0/25	25ml	B10	B10	10.64
FPC3/25	25ml	B14	B14	10.64
FPC3/50	50ml	B14	B14	10.90
FPC3/100	100ml	B14	B14	12.28

Pear Shaped Flasks, Claisen-Vigreux

Item	Capacity ml	Cone Size	Socket Size	Price Each
FPV4/25	25ml	B14	B14	13.42
FPV4/50	50ml	B14	B14	14.60
FPV4/100	100ml	B14	B14	14.60

Pear Shaped Flasks, Distillation

Item	Capacity ml	Cone Size	Socket Size	Price Each
FPD2/10	10ml	B14	B14	7.64
FPD2/25	25ml	B14	B14	6.96
FPD2/50	50ml	B14	B14	7.64
FPD2/100	100ml	B14	B14	8.80

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Flasks



Pear Shaped Flasks, Single Neck

Item	Capacity ml	Socket Size	Height mm	Dia mm	Price Each
FP5/B10	5ml	B10	65	25	4.14
FP5/B14	5ml	B14	70	25	3.88
FP10/B10	10ml	B10	80	30	4.32
FP10/B14	10ml	B14	85	30	4.12
FP25/B10	25ml	B10	90	40	4.14
FP25/B14	25ml	B14	100	40	4.74
FP50/B14	50ml	B14	115	50	3.42
FP50/B19	50ml	B19	130	50	4.74
FP100/B14	100ml	B14	135	60	5.26
FP100/B19	100ml	B19	140	60	5.84
FP100/B24	100ml	B24	150	60	7.38

Pear Shaped Flasks, Centre Neck & Angled Side Neck

Item	Capacity ml	Centre Socket	Side Socket	Height mm	Price Each
FP50/11	50ml	B14	B14	115	9.58

Pear Shaped Flasks, Centre Neck, Angled Side Neck & Parallel Side Neck

Item	Capacity ml	Centre Socket	Side Socket	Height mm	Price Each
FP25/000	25ml	B10	B10	115	15.76

Pear Shaped Flasks, Vigreux with Angled Neck

Item	Capacity ml	Main Socket	Side Socket	Side Cone	Price Each
FP7/10	10ml	B14	B10	B14	14.14
FP7/25	25ml	B14	B10	B14	14.14
FP7/50	50ml	B14	B10	B14	14.60

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Flasks



Rotary Evaporator Flasks

Item	Capacity ml	Socket Size	Price Each
RF50/B24	50ml	B24	9.46
RF50/B29	50ml	B29	9.46
RF100/B24	100ml	B24	10.88
RF100/B29	100ml	B29	10.88
RF250/B24	250ml	B24	11.98
RF250/B29	250ml	B29	11.98
RF500/B24	500ml	B24	12.48
RF500/B29	500ml	B29	12.48
RF1L/B24	1000ml	B24	13.08
RF1L/B29	1000ml	B29	14.66
RF2L/B29	2000ml	B29	20.86



Flasks, Round Bottom, Short Neck

Item	Capacity ml	Socket Size	Price Each
FRS5/B10	5ml	B10	2.33
FRS5/B14	5ml	B14	2.15
FRS10/B10	10ml	B10	2.71
FRS10/B14	10ml	B14	2.57
FRS25/B14	25ml	B14	2.36
FRS25/B19	25ml	B19	2.43
FRS25/B24	25ml	B24	2.63
FRS50/B14	50ml	B14	2.80
FRS50/B19	50ml	B19	2.83
FRS50/B24	50ml	B24	2.94
FRS50/B29	50ml	B29	4.01
FRS100/B14	100ml	B14	3.50
FRS100/B19	100ml	B19	3.43
FRS100/B24	100ml	B24	3.47
FRS100/B29	100ml	B29	3.80
FRS100/B34	100ml	B34	4.08
FRS150/B19	150ml	B19	3.78
FRS150/B24	150ml	B24	3.80
FRS150/B29	150ml	B29	4.15
FRS150/B34	150ml	B34	4.06
FRS250/B14	250ml	B14	4.17
FRS250/B19	250ml	B19	4.01
FRS250/B24	250ml	B24	3.82
FRS250/B29	250ml	B29	4.24
FRS250/B34	250ml	B34	4.34

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Flasks



Flasks, Round Bottom, Short Neck

Item	Capacity ml	Socket Size	Price Each
FRS500/B19	500ml	B19	4.41
FRS500/B24	500ml	B24	4.41
FRS500/B29	500ml	B29	5.11
FRS500/B34	500ml	B34	5.37
FRS500/B45	500ml	B45	7.04
FRS1L/B19	1000ml	B19	5.86
FRS1L/B24	1000ml	B24	4.88
FRS1L/B29	1000ml	B29	5.99
FRS1L/B34	1000ml	B34	6.69
FRS1L/B45	1000ml	B45	9.24
FRS1L/B55	1000ml	B55	9.68
FRS2L/B24	2000ml	B24	10.38
FRS2L/B29	2000ml	B29	9.94
FRS2L/B34	2000ml	B34	10.03
FRS3L/B24	3000ml	B24	12.32
FRS3L/B34	3000ml	B34	12.85
FRS5L/B29	5000ml	B29	19.37
FRS5L/B34	5000ml	B34	20.23
FRS5L/B45	5000ml	B45	19.04
FRS5L/B55	5000ml	B55	20.79



Flasks, Round Bottom, Two Necks, Angled Side Neck

Item	Capacity ml	Centre Neck Socket Size	Angled Side Neck Socket Size	Price Each
FR100/31A	100ml	B24	B14	6.95
FR100/32A	100ml	B24	B19	8.14
FR250/31A	250ml	B24	B14	7.26
FR250/32A	250ml	B24	B19	7.82
FR500/31A	500ml	B24	B14	7.98
FR500/32A	500ml	B24	B19	9.14
FR1L/31A	1000ml	B24	B14	9.56
FR1L/32A	1000ml	B24	B19	10.33
FR2L/52A	2000ml	B34	B19	15.66

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Flasks



Flasks, Round Bottom, Three Necks, Parallel Side Necks

Item	Capacity ml	Centre Neck Socket Size	Parallel Side Necks Socket Size	Price Each
FR250/322P	250ml	B24	B19	12.34
FR250/411P	250ml	B29	B14	11.18
FR500/322P	500ml	B24	B19	13.28
FR500/411P	500ml	B29	B14	14.63
FR500/422P	500ml	B29	B19	13.95
FR1L/322P	1000ml	B24	B19	14.67
FR1L/411P	1000ml	B29	B14	16.99
FR1L/522P	1000ml	B34	B19	15.84
FR1L/533P	1000ml	B34	B24	17.48
FR1L/722P	1000ml	B45	B19	16.89
FR2L/322P	2000ml	B24	B19	23.91
FR2L/422P	2000ml	B29	B19	19.60
FR2L/522P	2000ml	B34	B19	20.35
FR2L/533P	2000ml	B34	B24	21.88
FR3L/522P	3000ml	B34	B19	26.86
FR3L/533P	3000ml	B34	B24	27.88
FR5L/533P	5000ml	B34	B24	33.69



Flasks, Round Bottom, Three Necks, Angled Side Neck

Item	Capacity ml	Centre Neck Socket Size	Angled Side Necks Socket Size	Price Each
FR100/211A	100ml	B19	B14	10.38
FR100/311A	100ml	B24	B14	10.24
FR100/322A	100ml	B24	B19	11.66
FR250/222A	250ml	B19	B19	12.64
FR250/311A	250ml	B24	B14	10.24
FR250/322A	250ml	B24	B19	11.66
FR250/411A	250ml	B29	B14	10.82
FR500/311A	500ml	B24	B14	10.59
FR500/322A	500ml	B24	B19	12.46
FR500/411A	500ml	B29	B14	11.22
FR1L/311A	1000ml	B24	B14	12.22
FR1L/322A	1000ml	B24	B19	13.20
FR1L/522A	1000ml	B34	B19	15.14
FR2L/322A	2000ml	B24	B19	20.77
FR2L/522A	2000ml	B34	B19	19.96
FR3L/522A	3000ml	B34	B19	25.04

E&OE Prices Subject to Change Without Notice Feb 2012.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk

