

# Thermometers & Tubes

**SGL - Scientific Glass Laboratories Ltd** of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

*Providing Scientific Glassware to Laboratories for over 40 years*

**Scientific Glass  
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

[www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Tubes



## Test Tubes

Item	Capacity ml	Height mm	Dia mm	Socket Size	Price Each
TT6/0	2ml	65	14	B10	3.51
TT10/0	5ml	100	14	B10	3.73
TT10/1	11ml	100	18	B14	3.10
TT10/2	19ml	100	23	B19	3.87
TT12/1	15ml	125	18	B14	3.34
TT15/1	19ml	150	18	B14	2.99
TT15/2	31ml	150	23	B19	3.77
TT15/3	54ml	150	29	B24	3.61
TT20/3	75ml	200	29	B24	4.10
TT25/5	195ml	250	41	B34	9.46

## Tubes, Drying, Angled

Item	Cone Size	Price Each
DTS/0	B10	7.26
DTS/1	B14	7.26
DTS/2	B19	8.17
DTS/3	B24	9.66
DTS/4	B29	10.23

## Tubes / Air Inlets, Bent with Cone

Item	Stem Length	Cone Size	Price Each
MCS/1B	150mm	B14	6.22
MCS/2B	150mm	B19	6.86
MCS/3B	150mm	B24	7.78

## Tubes / Air Inlets, Straight with Cone

Item	Stem Length	Cone Size	Price Each
MCS/0	100mm	B10	5.33
MCS/1	100mm	B14	5.33
MCS/2	100mm	B19	5.58
MCS/3	100mm	B24	5.63
MCS/4	100mm	B29	6.61

E&OE Prices Subject to Change Without Notice March 2021.  
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Thermometer



## Thermometer Pockets

Item	Stem Length	Cone Size	Price Each
TPS/1	45mm	B14	5.56
TPS/2	200mm	B19	5.58

## Thermometers with Cone

Item	Temp Range	Cone Size	Price Each
T1/1	0 -110°C	B14	32.02
T1/2	0 -200°C	B14	32.02
T1/3	0 -300°C	B14	28.79
T1/4	0 -400°C	B14	30.23

E&OE Prices Subject to Change Without Notice March2021.  
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)



Scientific Glass Laboratories  
- Traditional Craftsmanship  
- Proven Quality

Product Literature available in this series:

- Adapters
- Beakers & Flasks
- Bends
- Chromatography Columns
- Cold Traps
- Condensers
- Dean & Stark
- Digestion Tubes
- Distillation
- Drechsel Bottles & Heads
- Extractors
- Flasks
- Funnels
- Funnels Dropping Separating
- Galileo Liquid Thermometer
- Pipettes, Flasks & Cylinders
- Screw Caps, Stoppers & Stopcocks
- Screwthreads, Discs & Joints
- Specimen Tubes & Vials
- Stirrers
- Tubes & Thermometer Pockets

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)

