

# Chromatography Columns

**SGL - Scientific Glass Laboratories Ltd** of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

*Providing Scientific Glassware to Laboratories for over 40 years*

**Scientific Glass  
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

[www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Chromatography Columns



## Chromatography Columns with Sintered Disc

Item	Sinter Dia	Socket Size	Length mm	Price Each
R1/10	10	14/23	100	32.75
R1/20	10	14/23	200	32.75
R1/30	10	14/23	300	32.75
R1/40	10	14/23	400	32.75
R2/10	20	19/26	100	33.76
R2/20	20	19/26	200	34.94
R2/30	20	19/26	300	34.94
R2/40	20	19/26	400	34.94
R2/50	20	19/26	500	34.94
R3/30	25	24/29	300	35.74
R3/40	25	24/29	400	35.74
R3/50	25	24/29	500	35.74
R5/40	30	34/35	400	41.20
R5/50	30	34/35	500	41.20
R6/40	40	40/45	400	55.89
R6/50	40	40/45	500	55.89

Columns 10mm to 25mm bore fitted with 2mm bore stopcocks.

Columns 30mm to 40mm bore fitted with 4mm bore stopcocks.

Glass or PTFE s/cs - please specify - /G or /P

Por 0, 1, 2, 3 and 4 available.

Add to code no. for example R1/10/0 etc.

## Chromatography Columns Reservoirs

Item	Size	Socket	Cone	Each
CCR-SRB24-100-CNB24	100ml	B24	B24	13.75
CCR-SRB24-250-CNB24	250ml	B24	B24	14.48
CCR-SRB24-500-CNB24	500ml	B24	B24	17.38
CCR-SRB24-1000-CNB24	1000ml	B24	B24	18.83
CCR-SRB24-2000-CNB24	2000ml	B24	B24	31.12

Other sizes & joints are available on request

## Chromatography Columns

A full range of bespoke columns are available on request.

Up to 130mm diameter tubing.

Columns without sinters are also available.

E&OE Prices Subject to Change Without Notice January 2023.  
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)

