

Adapters

SGL - Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

Providing Scientific Glassware to Laboratories for over 40 years

**Scientific Glass
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Adapters



Adapters, Cone / Rubber Tubing, right Angled Connection

Item	Cone Size	Price Each	With Thread Cap & Nozzle
ATA/0	B10	8.11	11.45
ATA/1	B14	7.52	9.45
ATA/2	B19	7.90	9.78
ATA/3	B24	7.97	9.83
ATA/4	B29	9.87	11.49
ATA/5	B34	11.41	12.83

For GL14 thread type add /T to code no.



Adapters, Cone / Rubber Tubing, Straight Connection

Item	Cone Size	Price Each
ATS/1	B14	7.52
ATS/2	B19	7.76
ATS/3	B24	7.97



Expansion Adapters

Item	Cone Size	Socket Size	Price Each
AX/10	B10	B14	9.21
AX/21	B14	B19	7.29
AX/31	B14	B24	6.33
AX/32	B19	B24	6.33
AX/41	B14	B29	9.78
AX/42	B19	B29	9.78
AX/43	B24	B29	8.05
AX/52	B19	B34	11.33
AX/53	B24	B34	9.69
AX/54	B29	B34	10.47
AX/64	B29	B40	14.29

E&OE Prices Subject to Change Without Notice January 2023.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Multiple, Two Necks Parallel

Item	Socket Size	Cone Size	Price Each
AMP/00	B10	B10	17.45
AMP/11	B14	B14	16.79
AMP/22	B19	B19	17.49
AMP/32	B19	B24	18.25
AMP/52	B19	B34	22.92
AMP/33	B24	B24	26.51
AMP/44	B29	B29	29.90

Adapters, Multiple, One Parallel & One 45 Degree Neck

Item	Socket Size	Cone Size	Price Each
AMA/22	B19	B19	17.15
AMA/32	B19	B24	17.15
AMA/33	B24	B24	17.15
AMA/43	B24	B29	17.93

Adapters, Multiple, Two Parallel & One 45 Degree Neck

Item	Socket Size	Cone Size	Price Each
AM/322	B19	B24	25.92
AM/422	B19	B29	25.92

Adapters, Receiver, Intermediate (Perkin) All Glass Taps

Item	Socket Size	Cone Size	Price Each
RIP/11	B14	B14	361.11
RIP/22	B19	B19	385.19
RIP/23	B19	B24	385.19
RIP/33	B24	B24	377.17

E&OE Prices Subject to Change Without Notice January 2023.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Vertical, Sloping

Item	Socket Size	Cone Size	Price Each
ARA/22B	B19	B19	16.70

Adapters, Vertical, Delivery

Item	Socket Size	Cone Size	Price Each
ARA/22	B19	B19	11.95

Adapters, Receiver, Vacuum (to take 4 x 100ml flasks)

Item	Socket Size	Cone Size	Price Each
ARA/21	B19	B14	41.83

Adapters, Receiver, Vertical with Vacuum Connection

Item	Socket Size	Cone Size	Price Each	With Thread Cap & Nozzle
RASV/22	B19	B19	18.84	22.18
RASV/33	B24	B24	18.84	22.18

For GL14 thread type add /T to code no.

E&OE Prices Subject to Change Without Notice January 2023.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Reduction Adapters

Item	Socket Size	Cone Size	Price Each
AR/01	B10	B14	5.65
AR/12	B14	B19	4.91
AR/13	B14	B24	5.38
AR/14	B14	B29	5.83
AR/23	B19	B24	5.12
AR/24	B19	B29	5.75
AR/25	B19	B34	6.13
AR/26	B19	B40	9.37
AR/27	B19	B45	10.38
AR/34	B24	B29	5.65
AR/35	B24	B34	5.91
AR/36	B24	B40	8.76
AR/37	B24	B45	9.19
AR/38	B24	B50	22.06
AR/39	B24	B55	24.38
AR/45	B29	B34	7.11
AR/47	B29	B45	10.88
AR/56	B34	B40	9.53
AR/58	B34	B50	22.06

Adapters, Cone / Screw Thread

Item	Cone Size	Thread Size	Dia Through Hole mm	Price Each
AST/1	B14	14	6.0 - 7.0	12.66
AST/2	B19	14	6.0 - 7.0	12.66
AST/3	B24	14	6.0 - 7.0	14.59
AST/4	B29	14	6.0 - 7.0	15.96
AST/5	B34	14	6.0 - 7.0	20.21
AST/6	B14	18	7.0 - 8.5	18.02
AST/7	B19	18	7.0 - 8.5	17.28
AST/8	B24	18	7.0 - 8.5	21.40
AST/9	B29	18	7.0 - 8.5	21.40
AST/11	B19	24	10.5 - 11.5	27.22
AST/12	B24	24	10.5 - 11.5	27.04
AST/13	B29	24	10.5 - 11.5	28.53
AST/14	B34	24	10.5 - 11.5	32.75

E&OE Prices Subject to Change Without Notice January 2023.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ

email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Cone / Socket with 'T' Connection

Item	Socket Size	Cone Size	Price Each	With Thread, Cap & Nozzle
ATR/11	B14	B14	9.18	12.52
ATR/12	B14	B19	9.49	12.83
ATR/13	B14	B24	10.22	13.57
ATR/22	B19	B19	12.67	16.02
ATR/23	B19	B24	12.67	16.02
ATR/24	B19	B29	10.68	14.03
ATR/25	B19	B34	11.03	14.39
ATR/33	B24	B24	12.71	16.05
ATR/44	B29	B29	13.48	16.79

For GL14 thread type add /T to code no.



Adapters, Socket / Rubber Tubing, Straight Connection

Item	Socket Size	Price Each	With Thread Cap & Nozzle
ATF/2	B19	9.49	12.83
ATF/3	B24	9.49	12.83
ATF/4	B29	9.78	12.83

For GL14 thread type add /T to code no.



Adapters, Swan Neck, Thread Size GL14

Item	Socket Size	Cone Size	Price Each
ASN/1	B14	B14	20.36
ASN/2	B19	B19	20.36
ASN/3	B19	B24	20.53
ASN/4	B24	B24	23.87

E&OE Prices Subject to Change Without Notice January 2023.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Adapters



Adapters, Stopcock, Cone / Rubber Tubing, Straight Connection

Item	Cone Size	Price Each	With Thread Cap & Nozzle
ASS/1	B14	27.31	30.66
ASS/2	B19	27.31	30.66
ASS/3	B24	27.59	30.95
ASS/4	B29	28.09	31.44

For GL14 thread type add /T to code no.
For P.T.F.E key add P to code no.

Adapters, Stopcock, Cone / Rubber Tubing, Bent Connection

Item	Cone Size	Price Each	With Thread Cap & Nozzle
ASA/1	B14	27.31	30.66
ASA/2	B19	27.31	30.66
ASA/3	B24	27.59	30.95
ASA/4	B29	28.09	31.44

For GL14 thread type add /T to code no.
For P.T.F.E key add P to code no.

Adapters, Stopcock, Socket / Rubber Tubing, Straight Connection

Item	Cone Size	Price Each	With Thread Cap & Nozzle
ASF/1	B14	27.31	30.66
ASF/2	B19	27.31	30.66
ASF/3	B24	28.09	31.44
ASF/2B	B19	27.59	30.95

For GL14 thread type add /T to code no.
For P.T.F.E key add P to code no.

Adapters, Straight Delivery

Item	Socket Size	Price Each
RAS/1	B14	7.61
RAS/2	B19	7.99
RAS/3	B24	8.40
RAS/4	B29	8.74
RAS/5	B34	9.06

E&OE Prices Subject to Change Without Notice January 2023.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk

