

# Condensers

**SGL - Scientific Glass Laboratories Ltd** of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- **Extensive Range of Laboratory Glassware**
- **Glassblown Bespoke Items**
- **Overnight Delivery Service**

*Providing Scientific Glassware to Laboratories for over 40 years*

**Scientific Glass  
Laboratories Ltd**

t. 01782 816237

f. 01782 575637

[www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Condensers



## Condensers, Air

Item	Surface Area cm <sup>2</sup>	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CA12/00	30	100	120	B10	B10	9.27
CA12/11	100	200	225	B14	B14	8.76
CA12/12	100	150	170	B19	B19	9.69
CA12/22	200	300	320	B19	B19	10.25
CA12/23	240	500	570	B24	B24	10.19
CA12/33	400	400	430	B19	B24	15.48
CA12/43	400	550	570	B24	B24	15.48
CA12/44	500	550	570	B29	B29	16.39

## Condensers, Allihn

Item	Surface Area cm <sup>2</sup>	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CAL20/2	130	200	285	-	B19	37.78
CAL20/3	130	200	285	-	B24	39.64
CAL20/4	130	200	295	-	B29	40.23
CAL20/5	130	200	310	-	B34	37.38
CAL20/6	130	200	310	-	B40	43.33
CAL20/8	150	200	310	-	B50	47.04
CAL30/33	180	30	400	B24	B24	50.29

For GL14 Thread Type add "T" to code no.



## Condensers, Coil (Graham)

Item	Surface Area cm <sup>2</sup>	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CC15/22	230	150	265	B19	B19	28.11
CC15/33	230	150	265	B24	B24	29.66
CC30/22	450	300	420	B19	B19	48.66
CC30/33	450	300	420	B24	B24	48.66
CC30/44	450	300	420	B29	B29	48.66
CCX15/4	230	150	260	-	B29	30.28
CCX15/5	230	150	260	-	B34	31.69
CCX20/6	350	200	300	-	B40	41.42
CCX20/7	400	200	315	-	B45	47.32
CCX20/8	400	200	315	-	B50	47.32

For GL14 Thread Type add /T to code no.



E&OE Prices Subject to Change Without Notice January 2023.  
Prices Shown GB Pounds Sterling. Ex VAT.

**Scientific Glass Laboratories Ltd**

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)



# Condensers



## Condensers, Double Surface (Davies)

Item	Surface Area cm <sup>2</sup>	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CDS15/11	180	150	240	B14	B14	36.27
CDS15/22	180	150	250	B19	B19	35.91
CDS15/33	230	150	250	B24	B24	40.25
CDS20/22	250	200	300	B19	B19	37.94
CDS20/33	300	200	300	B24	B24	40.71
CDS20/44	300	200	315	B29	B29	44.07
CDS20/3	300	200	290	-	B24	39.33

For GL14 Thread Type add "T" to code no.



## Condensers, Jacketed Coil

Item	Surface Area cm <sup>2</sup>	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CCJ15/22	430	150	320	B19	B19	70.58
CCJ15/33	430	150	320	B24	B24	70.58
CCJ15/44	430	150	320	B29	B29	70.60
CCJ15/05	430	150	300	-	B34	66.29

For GL14 Thread Type add "T" to code no.



## Condensers, Liebig

Item	Effective Length mm	Socket Size	Cone Size	Price Each	With Thread, Cap & Nozzle
CL60/00	60	B10	B10	21.64	28.34
CL15/11	150	B14	B14	15.46	22.15
CL20/22	200	B19	B19	16.73	23.43
CL20/33	200	B24	B24	20.01	26.71
CL40/22	400	B19	B19	23.38	30.08
CL40/33	400	B24	B24	23.38	30.16
CL40/44	400	B29	B29	26.14	32.84

For GL14 Thread Type add "T" to code no.

E&OE Prices Subject to Change Without Notice January 2023.  
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)



Scientific Glass Laboratories  
- Traditional Craftsmanship  
- Proven Quality

Product Literature available in this series:

- Adapters
- Beakers & Flasks
- Bends
- Chromatography Columns
- Cold Traps
- Condensers
- Dean & Stark
- Digestion Tubes
- Distillation
- Drechsel Bottles & Heads
- Extractors
- Flasks
- Funnels
- Funnels Dropping Separating
- Galileo Liquid Thermometer
- Pipettes, Flasks & Cylinders
- Screw Caps, Stoppers & Stopcocks
- Screwthreads, Discs & Joints
- Specimen Tubes & Vials
- Stirrers
- Tubes & Thermometer Pockets

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ  
email. [sales@scientificglass.co.uk](mailto:sales@scientificglass.co.uk) [www.scientificglass.co.uk](http://www.scientificglass.co.uk)

