

Condensers

SGL- Scientific Glass Laboratories Ltd of Stoke-on-Trent Staffordshire are one of the UK's largest specialist glass laboratories, providing a prompt delivery service to commercial and educational laboratories in the UK, Europe & Global Markets.

- Extensive Range of Laboratory Glassware
- Glassblown Bespoke Items
- Overnight Delivery Service

Providing Scientific Glassware to Laboratories for over 40 years

Scientific Glass
Laboratories Ltd

t. 01782 816237

f. 01782 575637

www.scientificglass.co.uk



Condensers



Condensers, Air

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CA12/00	30	100	120	B10	B10	9.73
CA12/11	100	200	225	B14	B14	9.20
CA12/12	100	150	170	B19	B14	10.17
CA12/22	200	300	320	B19	B19	10.76
CA12/23	240	500	570	B24	B24	10.70
CA12/33	400	400	430	B19	B24	16.25
CA12/43	400	550	570	B24	B24	16.25
CA12/44	500	550	570	B29	B29	17.21

Condensers, Allihn

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CAL20/2	130	200	285	-	B19	39.67
CAL20/3	130	200	285	-	B24	41.62
CAL20/4	130	200	295	-	B29	42.24
CAL20/5	130	200	310	-	B34	39.25
CAL20/6	130	200	310	-	B40	45.50
CAL20/8	150	200	310	-	B50	49.39
CAL30/33	180	30	400	B24	B24	52.80

For GL14 Thread Type add "T" to code no.

Condensers, Coil (Graham)

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CC15/22	230	150	265	B19	B19	29.52
CC15/33	230	150	265	B24	B24	31.14
CC30/22	450	300	420	B19	B19	51.09
CC30/33	450	300	420	B24	B24	51.09
CC30/44	450	300	420	B29	B29	51.09
CCX15/4	230	150	260	-	B29	31.79
CCX15/5	230	150	260	-	B34	33.27
CCX20/6	350	200	300	-	B40	43.49
CCX20/7	400	200	315	-	B45	49.69
CCX20/8	400	200	315	-	B50	49.69

For GL14 Thread Type add /T to code no.

E&OE Prices Subject to Change Without Notice January 2026.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ

email. sales@scientificglass.co.uk www.scientificglass.co.uk



Condensers



Condensers, Double Surface (Davies)

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CDS15/11	180	150	240	B14	B14	38.08
CDS15/22	180	150	250	B19	B19	37.71
CDS15/33	230	150	250	B24	B24	42.26
CDS20/22	250	200	300	B19	B19	39.84
CDS20/33	300	200	300	B24	B24	42.75
CDS20/44	300	200	315	B29	B29	46.27
CDS20/3	300	200	290	-	B24	41.30

For GL14 Thread Type add "T" to code no.



Condensers, Jacketed Coil

Item	Surface Area cm ²	Effective Length mm	Overall Length mm	Socket Size	Cone Size	Price Each
CCJ15/22	430	150	320	B19	B19	74.11
CCJ15/33	430	150	320	B24	B24	74.11
CCJ15/44	430	150	320	B29	B29	74.13
CCJ15/05	430	150	300	-	B34	69.60

For GL14 Thread Type add "T" to code no.



Condensers, Liebig

Item	Effective Length mm	Socket Size	Cone Size	Price Each	With Thread, Cap & Nozzle
CL60/00	60	B10	B10	22.72	29.76
CL15/11	150	B14	B14	16.23	23.26
CL20/22	200	B19	B19	17.57	24.60
CL20/33	200	B24	B24	21.01	28.05
CL40/22	400	B19	B19	24.55	31.58
CL40/33	400	B24	B24	24.55	31.67
CL40/44	400	B29	B29	27.45	34.48

For GL14 Thread Type add "T" to code no.

E&OE Prices Subject to Change Without Notice January 2026.
Prices Shown GB Pounds Sterling. Ex VAT.

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ
email. sales@scientificglass.co.uk www.scientificglass.co.uk



Scientific Glass Laboratories

- Traditional Craftsmanship

- Proven Quality

Product Literature available in this series:

- Adapters
- Beakers & Flasks
- Bends
- Chromatography Columns
- Cold Traps
- Condensers
- Dean & Stark
- Digestion Tubes
- Distillation
- Drechsel Bottles & Heads
- Extractors
- Flasks
- Funnels
- Funnels Dropping Separating
- Galileo Liquid Thermometer
- Pipettes, Flasks & Cylinders
- Screw Caps, Stoppers & Stopcocks
- Screwthreads, Discs & Joints
- Specimen Tubes & Vials
- Stirrers
- Tubes & Thermometer Pockets

Scientific Glass Laboratories Ltd

t. 01782 816237 f. 01782 575637

Canal Lane, Tunstall, Stoke-on-Trent, Staffordshire ST6 4PQ

email. sales@scientificglass.co.uk www.scientificglass.co.uk

